

()446-913

314

(031)260-9600 <http://www.gcri.co.kr>

(hklee@mail.gcri.co.kr)

: 08 - 01

2008 2 27

:

:

:

IFN-

1.

2.

IFN-

1. IFN-

	Heparin Tube	3 - (Gray, Red, Purple cap Tube)
	5 mL	1) 3 () 1mL 2) tube blood coating 10 , 3) , 16 37 incubator incubation
		<u>37 Incubator가</u> : incubation 37 Incubator 16 ~ 24 <u>37 Incubator가</u> : : (12)
	,	: ~ (12 16 가)

* : 3 3